

WEEKLY DEFIBRILLATOR INSPECTION

School _____ Defibrillator Location _____

Ensure the cabinet alarms are working and replace the two 9 volt batteries once a year. Follow the directions posted on the inside of the cabinet behind the defibrillator. Check to see if the green ready light is blinking on the defibrillator. If the green light is not blinking, or if the defibrillator is “chirping”, call maintenance immediately. Examine the defibrillator to ensure the unit has not been tampered with or the supplies have not been removed. Send completed form to Maintenance the last Monday of every month, and retain a copy for your inspection file.

[illegible]